



Certificate of Analysis

Product Name	Human CD8 Antibody (M8A16) [APC], mAb, Mouse
Size/Cat. No.	25 tests V03003-25T
Formulation	PBS, pH 7.4, containing gelatin and 0.03% ProClin300.
Storage	The antibody should be stored at 4°C and protected from prolonged exposure to light. Do not freeze.
Purification	The antibody was purified by affinity chromatography and conjugated with APC under optimal conditions. The solution is free of unconjugated APC.
Lot No.	PA181006
Concentration	40 µg/mL
Date of Manufacturing	2018/08/01
Expiry Date	2019/08/01
Quality control	Flow cytometric analysis of antibody surface-stained cells.



Human peripheral blood lymphocytes were stained with or without Human CD8 Antibody (M8A16) [APC].

Certified by:  Date: _____
Antibody Production Director