


Certificate of Analysis

Product Name	Human CD8 Antibody (M8A16) [Biotin], mAb, Mouse
Size/Cat. No.	100 tests V02704-100T
Formulation	PBS, pH 7.4, containing gelatin and 0.03% ProClin300.
Storage	The antibody should be stored at 4°C and protected from prolonged exposure to light. Do not freeze.
Purification	The antibody was purified by affinity chromatography and conjugated with APC under optimal conditions. The solution is free of unconjugated APC.
Lot No.	PB181007
Concentration	60 µg/mL
Date of Manufacturing	2018/09/06
Expiry Date	2019/09/06
Quality control	Flow cytometric analysis of antibody surface-stained cells.



Human peripheral blood lymphocytes were stained with or without Human CD8 Antibody (M8A16) [Biotin], followed by PE conjugated Streptavidin (405204, Biolegend).

Certified by:  Date: _____
Antibody Production Director