

Human Recombinant CD27 Stable Cell Line
Cat. No. M00610

Version 04282015

I. INTRODUCTION

Catalog Number: M00610

Cell Line Name: GS-H2/CD27

Gene Synonyms: S152, S152. LPFS2, T14, TNFRSF7

Expressed Gene: Codon Optimized from NM_001242.4

Host Cell: GS-H2

Quantity: Two vials of frozen cells (1×10^6 per vial)

Stability: 20 passages

Application: *in vitro* functional assay

Freeze Medium: 95% complete growth medium, 5% DMSO

Complete Growth Medium: MEM, 10% FBS

Culture Medium: MEM, 10% FBS, 2 $\mu\text{g/ml}$ Puromycin, 200 $\mu\text{g/ml}$ Hygromycin B

Mycoplasma Status: Negative

Functional Performance: Using CD70 (the ligand of CD27), Signal / Background (S/B) > 3

Storage: Liquid nitrogen immediately upon receipt

II. BACKGROUND

CD27 is a member of the tumor necrosis factor receptor superfamily. It is currently of interest to immunologists as a co-stimulatory immune checkpoint molecule. It binds to ligand CD70, and plays a key role in regulating B-cell activation and immunoglobulin synthesis. This receptor transduces signals that lead to the activation of NF- κ B and MAPK8/JNK. Varlilumab is an antibody that binds to CD27 and is an experimental cancer treatment.

REPRESENTATIVE DATA

- Protein Expression Validation




Figure 1. Flow cytometry analysis of CD27 protein expression in GS-H2/CD27 cells. Red: GS-H2, Blue: GS-H2/CD27.

- Functional Validation by *in vitro* Assay




Figure 2. Functional evaluation of GS-H2/CD27 cells by measuring CD70 induced IL-8 production.

III. THAWING AND SUBCULTURING

Thawing Protocol

1. Remove the vial from liquid nitrogen tank and thaw cells quickly in a 37°C water-bath.
2. Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol and transfer the cells to a 15 ml centrifuge tube containing 9 ml of complete growth medium.
3. Pellet cells by centrifugation at 200 x g for 5 min, and remove the medium.
4. Resuspend the cells in complete growth medium.
5. Transfer the cell suspension to a 10 cm dish with 10 ml of complete growth medium.
6. Grow the cells in incubator with 37°C, 5 % CO₂.
7. Add antibiotic in the following day.

Sub-culturing Protocol

1. Centrifuge the cells at 200 x g for 5min, and remove the medium.
2. Resuspend the cells in culture medium and add the cells suspension to new culture dish.
3. Grow the cells in incubator with 37°C, 5 % CO₂.

Subcultivation Ratio: A subcultivation ratio of 1:4 to 1:8 is recommended

Medium Renewal: Every 2 to 3 days

IV. REFERENCES

1. Prasad KV, Ao Z, Yoon Y, Wu MX, Rizk M, Jacquot S, Schlossman SF (Jun 1997). "CD27, a member of the tumor necrosis factor receptor family, induces apoptosis and binds to Siva, a proapoptotic protein". Proceedings of the National Academy of Sciences of the United States of America. 94 (12): 6346–51.
2. Yamamoto H, Kishimoto T, Minamoto S (Nov 1998). "NF-kappaB activation in CD27 signaling: involvement of TNF receptor-associated factors in its signaling and identification of functional region of CD27". Journal of Immunology. 161 (9): 4753–9.
3. Akiba H, Nakano H, Nishinaka S, Shindo M, Kobata T, Atsuta M, Morimoto C, Ware CF, Malinin NL, Wallach D, Yagita H, Okumura K (May 1998). "CD27, a member of the tumor necrosis factor receptor superfamily, activates NF-kappaB and stress-activated protein kinase/c-Jun N-terminal kinase via TRAF2, TRAF5, and NF-kappaB-inducing kinase". The Journal of Biological Chemistry. 273 (21): 13353–8.

GenScript USA Inc,
860 Centennial Ave.
Piscataway, NJ 08854
Toll-Free: 1-877-436-7274
Tel: 1-732-885-9188, Fax: 1-732-210-0262
Email: product@genscript.com
Web: <http://www.genscript.com>
For Research Use Only.

Limited Use License Agreement

This is a legal agreement between you (Licensee) and GenScript USA Inc. governing use of GenScript's stable cell line products and protocols provided to licensee. By purchasing and using the stable cell line, the buyer agrees to comply with the following terms and conditions of this label license and recognizes and agrees to such restrictions:

- 1) The products are not transferable and will be used at the site where they were purchased. Transfer to another site owned by buyer will be permitted only upon written request by buyer followed by subsequent written approval by GenScript.
- 2) The purchaser cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party.
- 3) The products sold by GenScript are for laboratory and animal research purposes only. The products are not to be used on humans, for consumption, or for any unlawful uses.

GenScript USA Inc. will not assert against the buyer a claim of infringement of patents owned or controlled by GenScript USA Inc. and claiming this product based upon the manufacture, use or sale of a clinical diagnostic, therapeutic and vaccine, or prophylactic product developed in research by the buyer in which this product or its components has been employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on the use of this product for other purposes, contact Marketing Department, GenScript USA Inc., GenScript 860 Centennial Ave Piscataway, NJ 08854. Phone: 1-732-885-9188. Fax: 1-732-210-0262. Email: marketing@genscript.com.