

**CD28 Expressing Cell Line**  
**Cat. No. M00611**

**Version 04282015**

---

## **I. INTRODUCTION**

Catalog Number: M00611

Cell Line Name: GS-J1

Quantity: Two vials of frozen cells ( $1 \times 10^6$  per vial)

Stability: 20 passages

Application: in vitro functional assay

Freeze Medium: 95% complete growth medium, 5% DMSO

Culture Medium: RPMI 1640, 10% FBS

Mycoplasma Status: Negative

Functional Performance: For Ipilimumab, Signal / Background (S/B) > 3


Storage: Liquid nitrogen immediately upon receipt

## **II. BACKGROUND**

GS-J1 cells are an immortalized line of human T lymphocyte cells that are used to study acute T cell leukemia, T cell signaling, and the expression of various chemokine receptors susceptible to viral entry, particularly HIV. GS-J1 cells are useful in science because of their ability to produce interleukin 2.

### III. REPRESENTATIVE DATA

- Validation by *in vitro* Functional Assay



**Figure 1.** Functional evaluation of GS-J1 by cell-based anti-CTLA4 activity assay. The EC<sub>50</sub> curves studies with Ipilimumab were shown in 96 and 384-well plate formats, respectively. Human IgG1 was used as a negative control. GS-J1 cells were co-cultured with GS-C1/CD80 together with CTLA-4 fusion protein and Ipilimumab.

### IV. THAWING AND SUBCULTURING

#### Thawing Protocol

1. Remove the vial from liquid nitrogen tank and thaw cells quickly in a 37°C water-bath.
2. Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol and transfer the cells to a 15 ml centrifuge tube containing 9 ml of complete growth medium.
3. Pellet cells by centrifugation at 200 x g for 5 min, and remove the medium.
4. Resuspend the cells in complete growth medium.
5. Transfer the cell suspension to a 10 cm dish with 10 ml of complete growth medium.
6. Grow the cells in incubator with 37°C, 5 % CO<sub>2</sub>.
7. Add antibiotic in the following day.

#### Sub-culturing Protocol

1. Centrifuge the cells at 200 x g for 5min, and remove the medium.
2. Resuspend the cells in culture medium and add the cells suspension to new culture dish.
3. Grow the cells in incubator with 37°C, 5 % CO<sub>2</sub>.

Subcultivation density: 2.5 × 10<sup>5</sup> cells/ml~4 × 10<sup>5</sup> cells/ml

Medium Renewal: Every 2 to 3 days

**GenScript USA Inc,**  
860 Centennial Ave.  
Piscataway, NJ 08854  
Toll-Free: 1-877-436-7274  
Tel: 1-732-885-9188, Fax: 1-732-210-0262  
Email: [product@genscript.com](mailto:product@genscript.com)  
Web: <http://www.genscript.com>  
**For Research Use Only.**

## Limited Use License Agreement

This is a legal agreement between you (Licensee) and GenScript USA Inc. governing use of GenScript's stable cell line products and protocols provided to licensee. By purchasing and using the stable cell line, the buyer agrees to comply with the following terms and conditions of this label license and recognizes and agrees to such restrictions:

- 1) The products are not transferable and will be used at the site where they were purchased. Transfer to another site owned by buyer will be permitted only upon written request by buyer followed by subsequent written approval by GenScript.
- 2) The purchaser cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party.
- 3) The products sold by GenScript are for laboratory and animal research purposes only. The products are not to be used on humans, for consumption, or for any unlawful uses.

GenScript USA Inc. will not assert against the buyer a claim of infringement of patents owned or controlled by GenScript USA Inc. and claiming this product based upon the manufacture, use or sale of a clinical diagnostic, therapeutic and vaccine, or prophylactic product developed in research by the buyer in which this product or its components has been employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on the use of this product for other purposes, contact Marketing Department, GenScript USA Inc., 120 Centennial Avenue, Piscataway, New Jersey 08840, U.S.A. Phone: 1-732-885-9188. Fax: 1-732-210-0262. Email: [marketing@genscript.com](mailto:marketing@genscript.com).