

His 标签 ELISA 检测试剂盒

Cat. No. L00436 / Cat. No. L01010

该试剂盒含有 His 标签蛋白检测所需的所有试剂和缓冲液，His 标签 ELISA 检测试剂盒基于竞争性 ELISA 法，将 His 标签蛋白 (M.W.12.7 kD) 预先包被于微孔板，加入已知浓度的 His 标签蛋白标准品或样本，与预包被的 His 标签蛋白竞争结合 His 标签蛋白抗体，通过洗板、显色等步骤可以检测 His 标签蛋白的表达或半定量 His 标签蛋白。该试剂盒同时提供 His 标签蛋白储存液，可供客户添加回收检验以验证试剂盒的准确性。

产品应用


检测 His 标签蛋白

筛选最佳表达条件

用已知蛋白作为标准品进行半定量

产品特点

- ▶ **检测范围宽**: 1 ng/mL-729 ng/mL
- ▶ **灵敏度高**: 可检测 1 ng/mL His 标签蛋白 (如图示)
- ▶ **多功能性**: N/C 端 His 标签蛋白，来自真核或原核表达系统的 Hexa-His (6×His)、Penta-His 以及 Tetra-His 标签蛋白
- ▶ **便捷**: 含有检测所需的所有试剂，1.5 小时即可完成整个操作
- ▶ **兼容性强**: 适用于细菌、酵母、哺乳动物细胞裂解液和细胞培养物上清等绝大多数系统



产品名称	规格	产品编号
His Tag ELISA Detection Kit	1 Kit	L00436
His Tag ELISA Detection Kit (384-Well Plate)	384 Tests	L01010
His Tag Antibody Plate (black, 96 wells, non-removable)	5 plates	L00440B
His Tag Antibody Plate (white, 96 wells, non-removable)	5 plates	L00440W
His Tag Antibody Plate (clear, 8 wells × 12 strips)	5 plates	L00440C

其他标签检测产品

检测标签	产品名称	规格	产品编号
GST	GST tag ELISA Detection Kit	1 kit	L00411
	DYKDDDDK Tag Antibody Plate (black, 96 wells, non-removable)	5 plates	L00455B
Flag	DYKDDDDK Tag Antibody Plate (white, 96 wells, non-removable)	5 plates	L00455W
	DYKDDDDK Tag Antibody Plate (clear, 8 wells × 12 strips)	5 plates	L00455C