

His标签ELISA检测试剂盒

Cat. No. L00436 / Cat. No. L01010

该试剂盒含有His标签蛋白检测所需的所有试剂和缓冲液，His标签ELISA检测试剂盒基于竞争性ELISA法，将His标签蛋白 (M.W.12.7 kD) 预先包被于微孔板，加入已知浓度的His标签蛋白标准品或样本，与预包被的His标签蛋白竞争结合His标签蛋白抗体，通过洗板、显色等步骤可以检测His标签蛋白的表达或半定量His标签蛋白。该试剂盒同时提供His标签蛋白储存液，可供客户添加回收检验以验证试剂盒的准确性。

产品应用

检测His标签蛋白

筛选最佳表达条件

用已知蛋白作为标准品进行半定量

产品特点

● 灵敏度



检测范围宽

1 ng/ml-729 ng/ml



灵敏度高

可检测1 ng/ml His标签蛋白（如图示）



多功能性

N/C端His标签蛋白，来自真核或原核表达系统的
Hexa-His (6 × His)、Penta-His以及Tetra-His标签蛋白




便捷

含有检测所需的所有试剂，1.5小时即可完成整个操作



兼容性强

适用于细菌、酵母、哺乳动物细胞裂解液和细胞培养物
上清等绝大多数系统



产品列表

产品名称	规格	产品编号
His Tag ELISA Detection Kit	1 Kit	L00436
His Tag ELISA Detection Kit (384-Well Plate)	384 Tests	L01010
His Tag Antibody Plate (black, 96 wells, nonremovable)	5 plates	L00440b
His Tag Antibody Plate (white, 96 wells, nonremovable)	5 plates	L00440w
His Tag Antibody Plate (clear, 8×12 strips)	5 plates	L00440c

其他标签检测产品

检测标签	产品名称	规格	产品编号
GST	GST tag ELISA Detection Kit	1 kit	L00411
	DYKDDDDK Tag Antibody Plate (Black, 96-well)	5 plates	L00455B
Flag	DYKDDDDK Tag Antibody Plate (White, 96-well)	5 plates	L00455W
	DYKDDDDK Tag Antibody Plate (Clear, 8×12 strip)	5 plates	L00455C