# Genome engineering of *Corynebacterium glutamicum* as a chassis for bioconversion of C1 feedstocks

### Yu WANG\*

Tianjin Institute of Industrial Biotechnology, Chinese Academy of Sciences, No. 32 West 7th Avenue, Tianjin 300308, P. R. China Email: wang\_y@tib.cas.cn

The present biomanufacturing process is highly depended on sugar-based resources. Recently, abundant natural gas supplies and feasible conversion of methane to methanol have made methanol a promising C1 feedstock with low cost and good availability for biomanufacturing. *Corynebacterium glutamicum* is a well-known industrial workhorse for amino acid production. Herein, facile yet robust genome editing tools were developed and applied for engineering *C. glutamicum* to serve as a promising chassis for bioconversion of C1 feedstocks.

To achieve efficient genome engineering of C. glutamicum, we developed a complex CRISPR/Cas9 based toolbox. Cas9 and gRNA expression cassettes were reconstituted to combat Cas9 toxicity and facilitate effective termination of gRNA transcription. CRISPR/Cas9-mediated gene deletion/insertion/replacement was achieved efficiently. Furthermore, Cas9 variants recognizing various PAMs and nucleotide (cytidine and adenine) deaminases were integrated into the toolbox and facilitated precise base editing and multiple-locus editing. The developed genome editing tools were then adapted to work with our robot-based integrative laboratory automation system, establishing a powerful BioFoundry platform. With this platform, we rationally designed and experimentally engineered C. glutamicum to serve as a promising chassis for methanol bioconversion. Specifically, pentose and methanol utilization pathways were reconstructed to engineer C. glutamicum as a methanol-dependent synthetic methylotroph. Its cell growth relies on co-utilization of methanol and xylose, and most notably methanol is an indispensable carbon source. Due to the methanol-dependent characteristic, adaptive laboratory evolution was successfully applied to improving methanol utilization. The evolved mutant showed a 40-fold increase in cell growth and utilized <sup>13</sup>C-methanol and xylose with a high mole ratio of 7:1 and <sup>13</sup>C-labeling of intracellular metabolites up to 63%. Methanol-dependent amino acid production was also achieved, demonstrating the potential in bioconversion of methanol into valuable chemicals. Genome resequencing combined with reverse metabolic engineering revealed key genetic factors involved in methanol utilization in C. glutamicum.

## Yu Wang Associate Professor Tianjin Institute of Industrial Biotechnology, Chinese Academy of Sciences

#### Education:

PhD, 2012 – 2016, Biology, Shanghai Jiao Tong University, China MSc, 2009 – 2012, Microbiology, Shandong University, China BSc, 2005 – 2009, Microbial Technology, Shandong University, China

#### Professional Career:

2016 – 2018 Assistant Professor, Tianjin Institute of Industrial Biotechnology, CAS. 2018 – Present Associate Professor, Tianjin Institute of Industrial Biotechnology, CAS.

#### **Research Interests:**

Synthetic Biology	Metabolic Engineering
Genome Engineering	C1 Feedstock Bioconversion

#### Selected publications

- 1. **Wang, Y.**, et al., 2019. Expanding targeting scope, editing window, and base transition capability of base editing in *Corynebacterium glutamicum*. Biotechnol Bioeng DOI:10.1002/bit.27121.
- 2. **Wang, Y.**, et al., 2018. MACBETH: multiplex automated *Corynebacterium glutamicum* base editing method. Metab Eng 47, 200–210.
- 3. **Wang, Y.**, et al., 2018. A novel *Corynebacterium glutamicum* L-glutamate exporter. Appl Environ Microbiol 84, e02691-17.
- 4. Tuyishime, P.<sup>#</sup>, **Wang, Y.**<sup>#</sup>, et al., 2018. Engineering *Corynebacterium glutamicum* for methanol-dependent growth and glutamate production. Metab Eng 49, 220–231.
- 5. Fan, L.<sup>#</sup>, **Wang, Y.**<sup>#</sup>, et al., 2018. Engineering artificial fusion proteins for enhanced methanol bioconversion. ChemBioChem 19, 2465–2471.
- Liu, J.<sup>#</sup>, Wang, Y.<sup>#</sup>, et al., 2018. Mutations on peptidoglycan synthesis gene *ponA* improve electro-transformation efficiency of *Corynebacterium glutamicum* ATCC 13869. Appl Environ Microbiol 84, e02225-18.
- Wang, Y., et al., 2017. Switch of metabolic status: redirecting metabolic flux for acetoin production from glycerol by activating a silent glycerol catabolism pathway. Metab Eng 39, 90–101.
- 8. **Wang, Y.**, et al., 2015. Production of C3 platform chemicals from CO<sub>2</sub> by genetically engineered cyanobacteria. Green Chem 17, 3100–3110.
- 9. **Wang, Y.**, et al., 2014. Glycerol dehydrogenase plays a dual role in glycerol metabolism and 2,3-butanediol formation in *Klebsiella pneumoniae*. J Biol Chem 289, 6080–6090.

