


Positive

# MS Spectrum

Intensity



## Sample Information

User : Grace  
Date Acquired : 07/07/2009 21:21:23  
Sample Name : HIV gag199-207  
Cat No. : RP20251  
Lot No. : CD10907  
Injection Volum : 1.000µL  
Theoretical MW : 1064.33  
Observed MW : 1063.80

Interface : ESI  
Nebulizing Gas Flow : 1.50 L/min  
CDL Temp : 250°C  
CDL Volt : -25.0V  
Block Temp : 200°C

Prerod Bias : +4.5kv  
Detector : 1.5kv  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH