


Positive

Intensity

# MS Spectrum



## Sample Information

User : Grace  
Date Acquired : 08/14/2009 9:09:57  
Sample Name : BMf-BH3  
Cat No. : RP20201  
Lot No. : CA05908  
Injection Volum : 1.000µL  
Theoretical MW : 3012.50  
Observed MW : 3013.05

Interface : ESI  
Nebulizing Gas Flow : 1.50 L/min  
CDL Temp : 250°C  
CDL Volt : -25.0V  
Block Temp : 200°C

Prerod Bias : +4.5kv  
Detector : 1.5kv  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH