

Sample Name : TAT peptide
 Cat. No. : RP20256
 Lot. No. : CB06908


Pump A : 0.065% trifluoroacetic in 100% water
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile
 Total Flow : 1 ml / min
 Wavelength : 220 nm

Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
35.00	Pumps	Pump A B.Conc	95	
35.01	Pumps	Pump A B.Conc	5	
45.00	Pumps	Pump A B.Conc	5	
45.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : AlltimaTM C18 4.6 x 250 mm



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	9.865	3761450	535113	97.388
2	10.144	32793	4445	0.849
3	10.450	4451	641	0.115
4	14.290	5319	1033	0.138
5	14.723	9314	1668	0.241
6	19.722	10721	1018	0.278
7	20.994	28925	4123	0.749
8	23.018	9359	1145	0.242
Total		3862330	549185	100.000