


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 08/23/2009 10:20:33
Sample Name : Protease - Activated Receptor - 2, amide
Cat No. : RP20334
Lot No. : CA14908
Injection Volum : 1.000μL
Theoretical MW : 614.80
Observed MW : 614.25

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH