


Positive

# MS Spectrum

Intensity



## Sample Information

User : Grace  
Date Acquired : 08/23/2009 10:21:38  
Sample Name : MPG, HIV related  
Cat No. : RP20335  
Lot No. : CB08908  
Injection Volum : 1.000µL  
Theoretical MW : 2766.20  
Observed MW : 2765.70

Interface : ESI  
Nebulizing Gas Flow : 1.50 L/min  
CDL Temp : 250°C  
CDL Volt : -25.0V  
Block Temp : 200°C

Prerod Bias : +4.5kv  
Detector : 1.5kv  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH