


MS Spectrum Graph



Date and Time :2009-9-5 9:55:02
Instrument No. :0200013
User :WN
Product Name :73769-1 KY37-NH2
Lot No. :73769001081709QX
Inj. Volume :1
MW :3905.36

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp: 250°C
Block Temp: 200°C
Probe bias: +4.5kv
Detector: 1.5kv
T.Flow: 0.2ml/min
B.conc: 50%H2O50%ACN