


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 09/19/2009 13:17:31
Sample Name : Kemptide
Cat No. : RP10287
Lot No. : CB05909
Injection Volum : 1.000μL
Theoretical MW : 771.90
Observed MW : 771.45

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH