


Positive

MS Spectrum

Intensity



Sample Information		Interface		Prerod Bias	:+4.5kv
User	:Grace	Nebulizing Gas Flow	:1.50 L/min	Detector	:1.5kv
Date Acquired	: 09/16/2009 12:30:55	CDL Temp	: 250°C	T.Flow	: 0.2ml/min
Sample Name	:HIV - 1 Rev (34 - 50)	CDL Volt	: -25.0V	B.conc	: 50%H2O/50%MeOH
Cat No.	: RP20362	Block Temp	: 200°C		
Lot No.	: CA10909				
Injection Volum	: 1.000µL				
Theoretical MW	:2437.75				
Observed MW	:2436.90				