


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 09/25/2009 10:55:51
Sample Name : hgp100(25-33)
Cat No. : RP20344
Lot No. : CC08909
Injection Volum : 1.000μL
Theoretical MW : 1155.30
Observed MW : 1155.20

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH