


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 01/06/2010 15:46:09
Sample Name : LL37 (Human)
Cat No. : RP13323
Lot No. : 85691003121409BF
Injection Volum : 1.000μL
Theoretical MW : 4493.28
Observed MW : 4493.60

Interface : ESI
Nebulizing Gas Flow : 4.21 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH