

Positive

MS Spectrum

Intensity

60


50

40

30

20

10



m/z

Sample Information

User : Grace
Date Acquired : 10/28/2009 13:47:00
Sample Name : Tyrosine Protein Kinase JAK 2 (pTyr8, 9)
Cat No. : RP20119
Lot No. : CE05910
Injection Volum : 1.000µL
Theoretical MW : 2082.17
Observed MW : 2082.30

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH