


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 11/26/2009 15:59:02
Sample Name : InsB9-23
Cat No. : RP20391
Lot No. : 84280002112509ZW
Injection Volum : 1.000µL
Theoretical MW : 1645.90
Observed MW : 1646.20

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH