


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 11/26/2009 15:59:59
Sample Name : OVA G4 peptide
Cat No. : RP20398
Lot No. : 84280008112509ZW
Injection Volum : 1.000µL
Theoretical MW : 906.08
Observed MW : 905.70

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH