


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 01/28/2010 12:18:13
Sample Name : Brain Natriuretic Peptide (BNP) (1-32), human
Cat No. : RP11119
Lot No. : 85691002121409ZWJ
Injection Volum : 1.000μL
Theoretical MW : 3464.10
Observed MW : 3463.50

Interface : ESI
Nebulizing Gas Flow : 4.21 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH