


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 01/22/2010 11:45:13
Sample Name : Competence-Stimulating Peptide-12261
Cat No. : RP20427
Lot No. : 88894002012910BF
Injection Volum : 1.000µL
Theoretical MW : 2153.50
Observed MW : 2152.80

Interface : ESI
Nebulizing Gas Flow : 4.21 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH