


Positive

MS Spectrum



Sample Information		Interface		Prerod Bias	
User	:Grace	Interface	:ESI	Prerod Bias	:+4.5kv
Date Acquired	: 05/10/2009 13:53:50	Nebulizing Gas Flow	:1.50 L/min	Detector	:1.5kv
Sample Name	:Kallikrein Inhibitor	CDL Temp	: 250°C	T.Flow	: 0.2ml/min
Cat No.	: RP10161	CDL Volt	: -25.0V	B.conc	: 50%H2O/50%MeOH
Lot No.	: CC02905	Block Temp	: 200°C		
Injection Volum	: 1.000µL				
Theoretical MW	: 773.90				
Observed MW	: 773.25				