


MS Spectrum

Positive



Sample Information Item

Date and Time : 2005-3-9 19:01:45
 User : FJS
 Sample : RP10516(DM-7)
 Inj. Volume : 1.000
 MW : 766.92

Cat. No.: RP10516

Lot No.: PX00511503

Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp: 250°C
 Block Temp: 200°C

Probe bias: +4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. conc: 50%H2O/50%ACN