


MS Spectrum

Positive



Sample Information Item

Date and Time : 2005-4-1 17:36:36
 User : xie
 Sample : RP10818 FF-4-NH2
 Inj. Volume : 1.000
 MW : 598.77

Cat. No.: RP10818

Lot No.: PJ00115506

Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp: 250℃
 Block Temp: 200℃

Probe bias: +4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. conc: 50%H2O/50%ACN