


MS Spectrum

Positive



Sample Information Item

Date and Time : 2005-3-9 18:46:44
 User : FJS
 Sample : RP10887(YM-5-NH2)
 Inj. Volume : 1.000
 MW : 572.69

Cat. No.: RP10887

Lot No.: PX00211503

Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp: 250°C
 Block Temp: 200°C

Probe bias: +4.5kv
 Detector: 1.5kv
 T Flow: 0.2ml/min
 B-conc: 50%H2O/50%ACN