


Positive

MS Spectrum



Sample Information		Interface		Prerod Bias	
User	:Grace	Nebulizing Gas Flow	:ESI	Detector	:+4.5kv
Date Acquired	: 03/02/2009 19:10:13	CDL Temp	:1.50 L/min	T.Flow	:1.5kv
Sample Name	:Bradykinin	CDL Volt	: 250°C	B.conc	: 0.2ml/min
Cat No.	: RP11129	Block Temp	: -25.0V		: 50%H2O/50%MeOH
Lot No.	: CC02903		: 200°C		
Injection Volum	: 1.000µL				
Theoretical MW	: 1060.20				
Observed MW	: 1060.60				