


# MS Spectrum

Positive



## Sample Information Item

Date and Time : 2005-4-26 12:03:37  
 User : FJS  
 Sample : RP11253(VF-6)  
 Inj. Volume : 1.000  
 MW : 774.93

Cat. No.: RP11253

Lot No.: PL00512507

Probe: ESI  
 Nebulizer Gas Flow: 1.5L/min  
 CDL: -20.0v  
 CDL Temp: 250°C  
 Block Temp: 200°C

Probe bias: +4.5kv  
 Detector: 1.5kv  
 T. Flow: 0.2ml/min  
 B. conc: 50%H<sub>2</sub>O/50%ACN