


Positive

MS Spectrum

Intensity



Sample Information

User : Grace
Date Acquired : 11/20/2009 9:21:12
Sample Name : Dermorphin Analog
Cat No. : RP13522
Lot No. : 83177004111309WJ
Injection Volum : 1.000μL
Theoretical MW : 901.43
Observed MW : 901.70

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH