


Summary:

Column : 250*4.6mm, Hypersil ODS-2
Solvent A : 0.1% TFA in 100% water
Solvent B : 0.1% TFA in 100% acetonitrile
Gradient :
 0.1min 90 10
 25 min 90 10
 25.01min 0 100
 30 min 0 100
Flow rate : 1.0ml/min

Cat. No.: RP10112

Lot No.: PL00208503



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		5.563	17683.477	151658.438	0.6437
2		6.082	840425.000	23409160.000	99.3563
TOTAL			858108.477	23560818.438	100.0000