

Sample Description:

Summary:

Product Name: RP10142(KR-13)
Column : 4.6×250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water

Cat. No.: RP10142


Lot No.: PJ00508506

Gradient

	A	B
0.01min	13%	87%
25.0min	25%	75%
25.1min	100%	0%
30.0min		STOP

Flow rate : 1.0ml/min

Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		7.975	1089.769	8506.801	0.3573
2		14.657	346.579	3416.401	0.1435
3		15.072	342.924	4174.155	0.1753
4		15.370	348.532	2495.752	0.1048
5		15.667	1647.266	22298.998	0.9366
6		16.193	179236.797	2335556.000	98.0986
7		16.193	554.068	4377.194	0.1839
Total			183565.935	2380825.302	100.0000