

Sample Description:
Summary:

Product Name: RP10398 (AY-36-1)
Column : 4.6×250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water


Gradient

	A	B
0.01min	28%	72%
25.0min	42%	58%
25.1min	100%	0%
30.0min		

STOP

Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP10398
Lot No.: PJ00501506



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		14.743	1532.066	10054.598	0.0641
2		15.665	2448.176	21002.900	0.1340
3		17.523	2987.991	31900.145	0.2035
4		19.287	2926.847	35009.051	0.2233
5		21.947	34246.246	374012.313	2.3857
6		22.285	12710.100	107667.125	0.6868
7		22.830	524287.844	14929674.000	95.2327
8		23.792	7339.568	167733.484	1.0699
Total			588478.838	15677053.615	100.0000