

Sample Description:

Product Name : RP10442(RL-6)  
 Column : 4.6\*250mm, Kromasil C18 5um  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient


Cat. No.: RP10442

Lot No.: PJ00806506

	A	B
0.01min	20%	80%
25.0min	30%	70%
25.1min	100%	0%
30.0min		

STOP

Flow rate : 1.0ml/min  
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		11.132	372.545	5589.550	0.0537
2		12.317	9001.148	116295.195	1.1179
3		12.982	3993.438	33239.535	0.3195
4		13.225	11191.767	218321.188	2.0986
5		13.633	320582.125	10020999.000	96.3270
6		15.843	189.845	1967.401	0.0189
7		20.780	418.792	6693.750	0.0643
Total			345749.660	10403105.619	100.0000