


Sample Description:
Summary:

Product Name : RP10474 (PEM-11-2)
Column : 4.6*250mm, Kromasil C18 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water

Cat. No.: RP10474
Lot No.: PJ00313505

Gradient
0.01min A 20% B 80%
25.0min 40% 60%
25.1min 100% 0%
30.0min STOP

Flow rate : 1.0ml/min
Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		13.365	1402.740	16506.545	0.4999
2		16.705	4825.549	68899.422	2.0865
3		17.102	1791.292	46303.535	1.4022
4		17.862	216802.500	3156500.250	95.5897
5		17.862	1745.302	13925.531	0.4217
Total			226567.384	3302135.283	100.0000