


Summary:
 Column : 250*4.6mm, Hypersil ODS-2
 Solvent A : 0.1% TFA in 100% water
 Solvent B : 0.1% TFA in 100% acetonitrile
 Gradient : A B
 0.1min 75 25
 25 min 50 50
 25.01min 0 100
 30 min 0 100
 Flow rate : 1.0ml/min

Cat. No.: RP10509
 Lot No.: PL00412507



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		5.280	490.466	3449.850	0.0367
2		6.472	1054.397	8650.000	0.0919
3		9.605	459.757	4496.302	0.0478
4		11.660	444509.375	9364237.000	99.5332
5		12.512	1628.110	27319.385	0.2904
TOTAL			448142.106	9408152.536	100.0000