


Summary:  
 Column : 250\*4.6mm, Hypersil ODS-2  
 Solvent A : 0.1% TFA in 100% water  
 Solvent B : 0.1% TFA in 100% acetonitrile  
 Gradient :  
           A          B  
   0.1min  75      25  
   25 min  50      50  
   25.01min 0     100  
   30 min   0     100  
 Flow rate : 1.0ml/min

Cat. No.: RP10606  
 Lot No.: PX00202502



# RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		10.162	441.692	4493.094	0.4695
2		10.773	68771.508	930469.000	97.2178
3		12.105	136.743	1593.099	0.1665
4		13.350	1424.478	20541.646	2.1462
TOTAL			70774.422	957096.839	100.0000