


Sample Description:  
Summary:

Product Name : RP10673 (PEG-10-NH2)  
Column : 4.6\*250mm, Kromasil C18 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient  
0.01min 25% 75%  
25.0min 40% 60%  
25.1min 100% 0%  
30.0min  
Flow rate : 1.0ml/min  
Wavelength : 220nm

Cat. No.: RP10673  
Lot No.: PJ00210505



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		14.770	4212.921	81852.398	1.1997
2		15.437	5015.334	87629.352	1.2844
3		16.537	242187.031	6653335.500	97.5160
Total			251415.286	6822817.250	100.0000