


Summary:  
 Column : 250\*4.6mm, Hypersil ODS-2  
 Solvent A : 0.1% TFA in 100% water  
 Solvent B : 0.1% TFA in 100% acetonitrile  
 Gradient :       A       B  
           0.1min 75       25  
           25 min 50       50  
           25.01min 0       100  
           30 min 0       100  
 Flow rate : 1.0ml/min

Cat. No.: RP10741  
 Lot No.: PX00102502



### RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		13.988	2080.768	26898.482	0.1259
2		14.632	16493.145	184690.953	0.8643
3		14.902	15402.757	208025.156	0.9735
4		15.390	982408.063	20707748.000	96.9067
5		15.390	15819.589	186762.422	0.8740
6		17.202	2888.792	54618.906	0.2556
TOTAL			1035093.112	21368743.920	100.0000