

Sample Description:


Product Name : SA-41-NH2-2  
 Column : 4.6\*250mm, Kromasil C18 5um  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient

Cat. No.: RP10983  
 Lot No.: PJ00306506

	A	B
0.01min	10%	68%
25.0min	53%	47%
25.1min	100%	0%

STOP

Flow rate : 1.0ml/min  
 Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		12.190	1953.191	26414.396	0.2184
2		13.847	2571.019	33902.613	0.2804
3		14.340	3372.183	37338.301	0.3088
4		15.732	478923.344	11741114.000	97.0985
5		15.732	17924.533	191532.234	1.5840
6		17.015	4599.745	61659.316	0.5099
Total			509344.015	12091960.861	100.0000