


Sample Description:

Summary:

Product Name: NF-32(0)
Column : 4.6×250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient
0.01min 20% 80%
25.0min 45% 55%
25.1min 100% 0%
30.0min STOP
Flow rate : 1.0ml/min
Wavelength : 220nm

Cat. No.: RP11121

Lot No.: PJ01425504



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		14.275	1870.339	16372.649	0.2324
2		16.598	327880.844	7028685.000	99.7676
Total			329751.183	7045057.649	100.0000