


Summary:  
 Column : 250\*4.6mm, Hypersil ODS-2  
 Solvent A : 0.1% TFA in 100% water  
 Solvent B : 0.1% TFA in 100% acetonitrile  
 Gradient :  
           A          B  
 0.1min 75          25  
 25 min 50          50  
 25.01min 0         100  
 30 min 0          100  
 Flow rate : 1.0ml/min

Cat. No.: RP11253  
 Lot No.: PL00512507



# RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		8.303	785.212	8968.781	0.1335
2		8.488	1309.466	11777.741	0.1753
3		9.268	5486.713	49229.352	0.7329
4		9.588	440413.875	6643366.000	98.8983
5		11.025	439.000	4032.296	0.0600
TOTAL			448434.266	6717374.170	100.0000