

## Sample Description:

## Summary:

Product Name : RP12784 (YR-11-NH2)  
Column : 4.6\*250mm, Hypersil-ODS2 5um  
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
Solvent B : 0.1% trifluoroacetic in 100% water  
Gradient

Cat. No.: RP12784


Lot No.: PJ00425504

	A	B
0.01min	7%	93%
25.0min	20%	80%
25.1min	100%	0%
30.0min		

STOP

Flow rate : 1.0ml/min

Wavelength : 220nm



## Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		17.240	5245.434	67266.727	0.6656
2		17.675	8773.724	169077.234	1.6731
3		18.332	237330.594	9869156.000	97.6612
Total			251349.751	10105499.961	100.0000