

Sample Description:


Product Name : RP12950 (BOC-FF-5)
Column : 4.6*250mm, Sinochrom ODS-BP 5um
Solvent A : 0.1% trifluoroacetic in 100% acetonitrile
Solvent B : 0.1% trifluoroacetic in 100% water
Gradient

Cat. No.: RP12950

Lot No.: PX00113507

0.01min A B
 50% 50%
25.0min 100% 0%
30.0min STOP

Flow rate : 1.0ml/min
Wavelength : 220nm



Results

Peak No.	Peak ID	Ret Time	Height	Area	Conc.
1		16.367	1350.582	8496.298	0.2135
2		18.070	614.032	4325.000	0.1087
3		18.703	1761.247	28436.195	0.7147
4		19.067	2026.437	13521.096	0.3398
5		19.735	443431.406	3865496.250	97.1481
6		22.393	5155.366	50241.031	1.2627
7		24.572	1152.373	8457.851	0.2126
Total			455491.444	3978973.721	100.0000