


Summary:

Structure : FQ-17  
 Column : 4.6×250mm, Kromasil  
 Solvent A : 0.1% trifluoroacetic in 100% acetonitrile  
 Solvent B : 0.1% trifluoroacetic in 100% water  
 Gradient :

Cat. No. : RP13084  
 Lot No. : PY00125501

	A	B
0.01min	15%	85%
25min	30%	70%
25.1min	100%	0%
30min	STOP	

Flow rate : 1.0 mL/min  
 Wavelength : 220nm



RESULTS

Peak No.	Peak ID	Ret Time	Height	Area	Conc
1		16.482	881.706	8750.400	0.2295
2		17.203	8628.705	117029.586	3.0697
3		17.865	213867.234	3627131.500	95.1392
4		19.217	5499.910	59535.695	1.5616
TOTAL			228877.555	3812447.182	100.0000