


Positive

Intensity

MS Spectrum



Sample Information

User :xu
Date Acquired : 2008-4-23 10:00:27
Sample Name :Flu-M1
Cat No. : RP19978
Injection Volum : 1.000μL
Theoretical MW :966.20
Observed MW :965.75

Interface :ESI
Nebulizing Gas Flow :1.50 L/min
CDL Temp : 250℃
CDL Volt : -25.0V
Block Temp : 200℃

Prerod Bias :+4.5kv
Detector :1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH