


Positive

MS Spectrum



Sample Information		Interface		Prerod Bias	:	+4.5kv
User	:xu	Nebulizing Gas Flow	:	Detector	:	1.5kv
Date Acquired	: 2008-4-23 9:59:17	CDL Temp	:	T.Flow	:	0.2ml/min
Sample Name	:Tetanus toxin (830-843)	CDL Volt	:	B.conc	:	50%H2O/50%MeOH
Cat No.	: RP19980	Block Temp	:			
Injection Volum	: 1.000µL					
Theoretical MW	:1611.90					
Observed MW	:1612.20					