


Positive

MS Spectrum



Sample Information		Interface	:ESI	Prerod Bias	:+4.5kv
User	:Grace	Nebulizing Gas Flow	:1.50 L/min	Detector	:1.5kv
Date Acquired	: 02/15/2009 11:57:06	CDL Temp	: 250℃	T.Flow	: 0.2ml/min
Sample Name	:Beta-Amyloid(1-17)	CDL Volt	: -25.0V	B.conc	: 50%H2O/50%MeOH
Cat No.	: RP20161	Block Temp	: 200℃		
Lot No.	: CA06902				
Injection Volum	: 1.000μL				
Theoretical MW	: 2068.20				
Observed MW	: 2068.40				