



Sample Information		Interface		Prerod Bias	:+4.5kv
User	:Grace	Nebulizing Gas Flow	:1.50 L/min	Detector	:1.5kv
Date Acquired	: 03/12/2009 11:29:46	CDL Temp	: 250°C	T.Flow	: 0.2ml/min
Sample Name	:Glucagon-Like Peptide(GLP)II,human	CDL Volt	: -25.0V	B.conc	: 50%H2O/50%MeOH
Cat No.	: RP10769	Block Temp	: 200°C		
Lot No.	: CC04903				
Injection Volum	: 1.000µL				
Theoretical MW	: 3766.20				
Observed MW	: 3765.45				