


Positive

MS Spectrum

Intensity



Sample Information

User : Grace  
Date Acquired : 04/30/2009 11:17:09  
Sample Name : Beta-Amyloid(1-34)  
Cat No. : RP20140  
Lot No. : CD01905  
Injection Volum : 1.000µL  
Theoretical MW : 3787.20  
Observed MW : 3786.60

Interface : ESI  
Nebulizing Gas Flow : 1.50 L/min  
CDL Temp : 250°C  
CDL Volt : -25.0V  
Block Temp : 200°C

Prerod Bias : +4.5kv  
Detector : 1.5kv  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH