

# His-tag antibody [Biotin], pAb, Rabbit

Cat. No.: A00205

## Overview

<b>Specificity</b>	GenScript His-tag antibody [Biotin], pAb, Rabbit specifically reacts with fusion proteins containing His tagged fusion proteins.
<b>Host Species</b>	Rabbit
<b>Immunogen</b>	His epitope tag peptide HHHHHH conjugated to KLH.
<b>Conjugate</b>	Biotin

## Applications


Working concentrations for specific applications should be determined by the investigator. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

<b>Application</b>	<b>Recommended Usage</b>
Western Blot	0.2-1 µg/ml

## Properties

<b>Form</b>	Lyophilized
<b>Storage Buffer</b>	Lyophilized with PBS, pH 7.4, containing 1% BSA and 0.02% sodium azide
<b>Reconstitution</b>	Reconstitute the lyophilized powder with deionized water (or equivalent) to a final concentration of 0.5 mg/ml.
<b>Storage Instructions</b>	The lyophilized product remains stable for 1 year at -20°C from the date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8°C or for 12 months at -20°C or below. Avoid repeated freezing and thawing cycles.
<b>Purification</b>	Affinity chromatography
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	Not applicable

## Examples



Western blot analysis of His tagged fusion protein using His-tag antibody [Biotin], pAb, Rabbit (GenScript, A00205, 1 µg/ml). The signal was developed with Streptavidin-HRP (GenScript M00091).

Predicted Size: 27 kD

Observed Size: 27 kD

## Background

**Target Background :** GenScript His-tag antibody [Biotin], pAb, Rabbit is the biotin conjugated Rabbit anti His antibody. The antibody recognizes fusion protein with His tag which is located at either C-terminal, or internal regions of the protein. It is suitable for detecting the expression level of His tagged proteins. GenScript His-tag antibody [Biotin], pAb, Rabbit is the biotin conjugated Rabbit anti His antibody. It is suitable for detecting the expression level of His tagged proteins.

**Synonyms :** His-tag antibody [Biotin], pAb, Rabbit;

**For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.**

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China